

## UPGRADE KIT

### THERATRON® EQUINOX™ EXTERNAL BEAM THERAPY SYSTEM

Best® Theratronics, the world's leading supplier of Gamma Teletherapy equipment, is pleased to announce the immediate availability of Theratron® Equinox™ External Beam Therapy System UPGRADE KIT!

Resulting from numerous requests from our valued Customers, the Equinox™ Upgrade provides for existing Theratron® 780C, 780E, Elite 80 and Phoenix™ units to be upgraded to modern, state-of-the-art technology in a timely and cost-attractive way.


#### Upgrade Features

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws (optional)
- R&V System Ready (optional)
- Service Log Files
- On-Board Verification
- Motorized Wedge (optional)
- Collision Detection (optional)

#### Upgrade Includes

- Complete removal of old controls and electronics, and installation of new control system and covers
- Replacing the old collimator system with the new Equinox Collimator
- Replacing old treatment table with new Avanza™ Table
- Retaining head rotation capability





**From this...**



**...to this!**

This upgrade allows older units to be upgraded into current units and has the possibility of future improvements such as:

- MLC for Conformal RT (replacement of shielding blocks)
- MLC for "Sliding Window" IMRT
- OBI (On Board Imaging)

This upgrade comes with two different couch options including:

**Avanza™ Patient Positioning Table:** Fully integrated with the Theratron® Equinox™, the Avanza™ provides accurate and reproducible patient positioning, along with efficient and comfortable patient set-up.

**Table 27 i:** With its five basic motions, isocentric, rotation, vertical, lateral, longitudinal and top rotation, the 27 i is a patient support table that offers a lower-cost, more 'hands-on' option for patient positioning.

For more information on the new Theratron® Equinox™ External Beam Teletherapy Upgrade Kit, please contact:

Customer Service Team in Canada: **613-591-2100, ext.1080**

Email: [service@theratronics.ca](mailto:service@theratronics.ca)

Write to: **Best® Theratronics Customer Support, 413 March Road, Ottawa, Ontario, K2K 0E4 Canada**

BT/MF 8015 TEUK (1) v.10.2009

