

Gammacell® 3000 Elan Source Options & CDR



GC 3000 Elan	Total Nominal Source Activity (+/- 20%)		Central Dose Rate (+/- 20%)		Estimated Timer for 25 Gy Central	Estimated Timer for 33 Gy Central
	Model	Ci	TBq	Gy/min	rad/min	Min
C	990	36.7	3.4	340	7.4	9.7
I	1450	53.7	5.0	500	5.0	6.6
II	2900	107.4	9.75	975	2.6	3.4

413 March Road, Ottawa, ON K2K 0E4 Canada
 phone 613 591 2100 fax 613 591 6627 www.theratronics.ca