

Gammacell® 3000 Elan Source Options & CDR



Model	Total Nominal Source Activity (+/- 15%)		Central Dose Rate (+/- 15%)		Estimated Timer for 25 Gy Central	Estimated Timer for 25 Gy Minimum
	Ci	TBq	Gy/min	rad/min	minutes	minutes
C	900	33.3	3.10	310	8.06	10.65
I	1300	48.1	4.48	448	5.58	7.37
II	2600	96.2	8.74	874	2.86	3.78

413 March Road, Ottawa, ON K2K 0E4 Canada
 phone 613 591 2100 fax 613 591 6627 www.theratronics.ca