


## Rad Disk Ultra™ With Absolute Filtration

### Features:

- ▶ Specifically designed for use in any Gammacell® 40 irradiator
- ▶ All autoclavable, polycarbonate construction mitigates attenuation of source and allows germ-free rodent transport from cage to Gammacell® and return, using standard micro-isolation techniques
- ▶ Attenuation factor of 2% (verified by Fricke dosimetry)
- ▶ Rodent divider system allows for any combination of two through eight compartments
- ▶ Proprietary, triple-layered absolute filter with protective screening maintains maximum protection for the rodent or material contents of the Rad Disk Ultra™
- ▶ Gammacell® 40's only manufacturer-certified rodent holder



BT/MF 8016 GC40E RDU (3)

v.02.2011

