

MOSFET Portable & Wireless Dosimetry Systems for Quality Assurance

- In-vivo measurements
- Simplifies dosimetry — minimizes QA
- IMRT QA & phantom work
- Brachytherapy
- IGRT / Tomotherapy
- Radiology
- Intracavitary measurements



Dose obtained in real-time



Best Theratronics GammaBeam Systems have been saving lives of cancer patients for over 60 years around the world!



Best™ GammaBeam 500 Total Body Irradiator System

The Total Body Irradiator GammaBeam 500 is a teletherapy unit designed to produce a large fixed rectangular radiation field at an extended source-to-skin distance in order to deliver total body irradiation. The unit can also be used for research and dosimetry purposes. Features include: dose rate up to 40 cGy/min at 2.5 m in air, record & verify and imaging capabilities, graphical touch screen data entry interface and motorized vertical motion of the head.

TeamBest Global Companies © 2021

Best Medical Canada (MOSFET orders) tel: 613 591 2100 877 668 6636 www.mosfet.ca

Best Theratronics tel: 613 591 2100 866 792 8598 www.theratronics.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA