

## Best™ Integrated Radiotherapy Solutions

TEAMBEST GLOBAL®



www.teambest.com

### Best GammaBeam™ 100/300 Equinox Teletherapy System with IMRT & IGRT



GammaBeam™  
saving lives

With **NEW** Multi-Leaf  
Collimator for 80 and 100 cm  
SAD units — IMRT, IGRT, SRS, SBRT  
and Tomotherapy capable with ActiveRx



### Best Avanza™ 6D Patient Positioning Table

The new Avanza 6D table provides adjustment in all 6 dimensions, including pitch and roll. This option will greatly improve the patient set-up for all treatments.



Best HDR Remote Afterloader




Best Treatment Planning Systems for LDR, HDR Brachytherapy, SRS, SBRT, Tomotherapy and Teletherapy with ActiveRx for IMRT



Advanced  
Keyboard  
Technology

Best Sonalis™  
Ultrasound  
Imaging System



### **NEW!** Best Sub-Compact Model 200 Self-Shielded Cyclotron with Optional Second Chemistry Module (Patent Pending)

- Low energy compact system — can be placed next to PET/CT
- Integrated self-shielded cyclotron and chemistry module
- Easy to operate push-button graphic interface
- Capable of producing: <sup>18</sup>FDG, Na<sup>18</sup>F, <sup>18</sup>F-MISO, <sup>18</sup>FLT, <sup>18</sup>F-Choline, <sup>18</sup>F-DOPA, <sup>18</sup>F-PSMA, <sup>11</sup>C, <sup>13</sup>N, <sup>68</sup>Ga and more!
- Ideal for Nuclear Cardiology and other Medical Applications

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

\* Certain products shown are not available for sale currently.