


# Best<sup>TM</sup> ABT Molecular Imaging

## BG-75 BIOMARKER GENERATOR



The BG-75 Biomarker Generator integrates a compact mini-cyclotron, kit based micro-chemistry, and automated quality control, simplifying in-house production of  $^{18}\text{F}$ -FDG and advanced biomarkers.

- Push button graphic interface
- Kit based chemistry
- Single or batch dose production
- Final dose delivery to syringe or vial (option)
- Automated quality control testing
- Integrated cyclotron & chemistry self-shielding
- Complete production lab in a 30 x 30 meter area



TeamBest Companies © 2019

Best ABT Molecular Imaging 3024 Topside Business Park Drive, Louisville, TN 37777 USA  
tel: 865 982 0098 800 336 4970 [www.bestabt.com](http://www.bestabt.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

