

Best™ Integrated Brachytherapy Solutions

Why Medical Physicists, Radiation Oncologists and Urologists Prefer Best™

Go for the Best™ for your patients!


“A unique aspect of the Best™ seed is its 5 mm length. Best also provides 5 mm laser-cut spacers, which do not have flared edges. With the Best® seed/spacer combination, we can be more confident of accurate seed placement.”¹

“The source strengths of Best™ seeds are uniformly distributed and the values always fall within the $\pm 3\%$ criterion.”³

Best Iodine-125 Seed



Best Palladium-103 Seed



“Best™ Model 2301 is the most isotropic of all the I-125 sources currently available. The source is a double encapsulated source.”²

“Best™ Model 2301 seeds are of superior design as they show up more distinctly in transverse ultrasound and fluoroscopic images.”³

Best Fiducial Markers



“Best™ Model 2301 characteristic x-rays have no significant effect on the photon spectrum and the source emits an essentially pure I-125 decay scheme spectrum.”^{2†}

“Delivery of the specific Sk for any given date, delivered on the date requested prior to the procedure. Extremely high level of flexibility, e.g. seed date changes, activity changes, etc. Prompt, courteous customer service, and low price.”⁴

¹ Communication from J. Limmer, UW Cancer Center Wausau Hospital, Wausau, Wisconsin – Dated December 14, 2000. ² B.H. Heintz, R.E. Wallace & J.M. Hevezi, Med. Phys. 28(4) 671-682 (2001).
³ Communication from Khai Lai, Kreiton Oncology Services Shoreline, Washington – Dated March 8, 2001. ⁴ Jeff Limmer, Medical Physicist.

You don't go to a shoe store and cut your feet to fit the shoe they have in stock ... then why compromise when it comes to improving the clinical outcome for your patients?

Best™ Seeds have the best visualization in ultrasound, fluoro, x-ray, CT and MRI, facilitating real-time dosimetry!

TEAMBEST GLOBAL®



www.teambest.com

Best Medical is the *only* company that makes custom seeds and strands to your exact specifications— shipped within 24 hours, 7 days a week, sterile & non-sterile!

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com