

Best™ Integrated Complete Brachytherapy Solutions




Best™ HDR Remote Afterloader

BestMick Applicator

BestMick 10/15/20/25 Cartridges



Best™ Kobold Tandem & Ring Applicator



Best™ Kobold Fletcher Tandem & Ovoid Applicator



Best™ Kobold Henschke Tandem & Ovoid Applicator




Best™ Cyber Sonalis™ Ultrasound Imaging System



Best™ Treatment Planning System



Best™ Stepper/Stabilizer



Best™ Esophageal Brachytherapy Applicator



Best™ Breast Double-Balloon Brachytherapy Applicator

Best™ Flexi Needles



15G Flexi Needle



16G Flexi Needle




17G Flexi Needle



15G Flexi Needle with SS Sleeve (Click Fit)



Flexi Needle with Female Hub



Best™ X-ray Marker Cables


Best™ vascular
Novoste

Vascular Brachytherapy for In-stent Restenosis

VBT Designed for the Cath Lab



- Convenient & Minimal Dosimetry Calculations
- Portable & Reusable
- Compatible with 6F Guide & 7F Guide-Extending Catheters
- Short Treatment Time



Available in 30 mm, 40 mm & 60 mm Source Train Lengths

*Certain products shown are not available for sale currently.

Best™ medical international
A TEAMBEST GLOBAL COMPANY


Best Medical International, Inc. • 7643 Fullerton Road • Springfield, VA 22153 USA
tel: 703 451 2378 800 336 4970 www.teambest.com

TeamBest Global Companies ©2025

Best™ Prostate, Rectal, GYN and Stepper/Stabilizer Templates




Best Medical offers a wide variety of disposable templates for HDR and LDR applications. Specify hole size when ordering (15, 16 or 17 gauge). These products should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.




Modified Rectal Template




Modified GYN Template



GYN Template




LDR Prostate Template




HDR Prostate Template



GYN Template



Grid Template for B-K type Stepper/Stabilizer*



Rectal Template

*BMI provides grid templates for all makes and models of steppers/stabilizers. **All dimensions are approximate.

Best™ Dosimetry Services

Portable Dosimetry System



mobileMOSFET Wireless Dosimetry System


MOSFET Array for Brachytherapy Procedures



Best™ Syed Multi-Channel Vaginal Applicator

This novel HDR applicator consists of a disposable central insert and an outer sleeve along with a positioning collar. *Source channels travel fully to the distal end, curving around to give more dose to the apex area, allowing for optimum dose delivery.*

Additionally, the outer sleeve provides a 5 mm wall, to reduce dose to the vaginal mucosa, with a 2 mm wall at the distal end of the applicator to allow for greater dose delivery. The applicator is 18 cm in length with diameters ranging from 2.0 cm to 3.5 cm.**




Best™ Brachytherapy Kit




5 Implant Needles



5 Stylets with Hubs



5 Single Leader Catheters**



5 Half Moon Buttons*



5 Red Caps



5 Friction Cuffs

* Half Moon Buttons are available in either radiopaque (pictured above) or clear nylon.

**Catheter tubes are available in purple, green, yellow, clear or blue with either radiopaque or clear nylon buttons.