



## Best™ Integrated Brachytherapy Solutions

More than 60 years have gone into making Brachytherapy more precise, with very few side effects or complications, and excellent clinical outcomes in the treatment of localized cancers.


TeamBest manufactures and provides comprehensive medical devices and accessories needed for a fully integrated Brachytherapy program including: providing assistance setting up, licensing, procurement of equipment, proctoring, reimbursement consulting, patient education and marketing.

In addition, TeamBest can make Brachytherapy treatment more affordable by charging per procedure!

Our Brachytherapy products are manufactured with pride in the USA – providing everything you need from One World • One Source • The Best!

- Ultrasound Systems
- HDR Afterloaders/Accessories
- LDR/HDR Treatment Planning Systems
- Stepper/Stabilizers
- LDR/HDR Templates
- Radioactive Seeds & Markers
- Needles & Applicators...and more!

Best Medical International, Inc. © 2020



Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
tel: 703 451 2378 800 336 4970 [www.besttotalsolutions.com](http://www.besttotalsolutions.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

\* Certain products shown are not available for sale currently.

## Best™ Fiducial Marker Kits



Best Medical's Gold Fiducial Markers are available in a variety of diameters and lengths to meet your needs for soft tissue localization such as prostate, liver, pancreas, lung, and breast.

- Best™ Gold Fiducial Markers are available in tissue-absorbable strands
- Standard strand consists of two markers at 2 cm apart center to center
- Special marker strands with custom spacing are also available
- Excellent visualization under X-ray, Fluoro, CT, MRI, Ultrasound and other imaging modalities
- Sold individually, stranded, and in pre-loaded, pre-sterilized needle sets
- Can be used for Brachytherapy, IMRT and IGRT

# Best™ Brachytherapy Kit



5 Implant Needles



5 Stylets with Hubs




5 Friction Cuffs



5 Half Moon Buttons\*



5 Red Caps



5 Single Leader Catheters\*\*



\* Half Moon Buttons are available in radiopaque (pictured) or clear nylon.


\*\* Catheter tubes are available in 4 colors (clear, blue, green or purple) with either radiopaque or clear nylon buttons.

## Best™ Multi-Channel Vaginal Applicator




This novel HDR applicator consists of a disposable central insert and an outer sleeve along with a positioning collar. Source channels in the central insert travel fully to the distal end to allow for optimum dose delivery.

Additionally, the outer sleeve provides a 5 mm wall, to reduce dose to the vaginal mucosa, with a 2 mm wall at the distal end of the applicator to allow for greater dose delivery. The applicator is 18 cm in length with diameters ranging from 2.0 cm to 3.5 cm.




## Best™ Brachytherapy Templates


Best Medical offers a wide variety of disposable templates for HDR and LDR applications. Specify hole size when ordering (15, 16 or 17 gauge). These products should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.




GYN  
Template  
(Disposable)




MODIFIED  
GYN  
Template  
(Disposable)




GYN  
Template



Prostate  
Template  
(Disposable)




Rectal  
Template  
(Disposable)



MODIFIED  
Rectal  
Template  
(Disposable)


## Best™ X-ray Marker Cables

Best X-ray Marker Cables are indicated for use in HDR brachytherapy procedures. All catheters are available in custom lengths, with or without cable/cable dummy ribbon.




## Best™ Prostate Templates

Best Medical International offers a variety of prostate grid templates for all makes and models of steppers, stabilizers and ultrasound systems.



Grid Template  
(Disposable)  
with Holder  
Bracket



Grid Template for  
B-K type (Disposable)