


Best[™] GammaBeam[™] 300-100 CM Equinox[™] Teletherapy System with Avanza 6D Patient Positioning Table



NEW!

Multi-Leaf Collimator (built in or external) for 80 and 100 cm SAD units—IMRT, IGRT, SRS, SBRT and Tomotherapy capable w/ActiveRx

UPGRADE KIT

for all old Theratron units, 80 or 100 cm including IMRT capabilities w/built in or external MLC



UPGRADE includes:

- Removing all old controls, electronics and installing a new control system and covers
- Replacing the old collimator system with the new Equinox collimator
- Replacing the old treatment table with the new Avanza[™] Table
- Retaining the head rotation capability is optional

UPGRADE features:

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws
- R&V System Ready
- Service Log Files
- On-Board Verification
- Motorized Wedge
- Collision Detection