



NEW! Best Model 200 Sub-Compact Self-Shielded Cyclotron with Optional Second Chemistry Module

- Low energy compact system — can be placed next to PET/CT
- Integrated self-shielded cyclotron and chemistry module
- Easy to operate push-button graphic interface
- Capable of producing: ^{18}F FDG, Na^{18}F , ^{18}F -MISO, ^{18}F FLT, ^{18}F -Choline, ^{18}F -DOPA, ^{18}F -PSMA, ^{11}C , ^{13}N , ^{68}Ga and more!
- Ideal for Nuclear Cardiology/ Oncology and other Medical Applications



NEW! Introducing a Novel Target for Best Model 200 Sub-Compact Self-Shielded Cyclotron (Patent Pending)

- Improved and easier clamp design
- Improved cooling
- Capable of double the beam current as original design
- Increased yield of Radiopharmaceuticals per run
- New material for better heat transfer

TeamBest Global Companies © 2022

Best Cyclotron Systems 3024 Topside Business Park Drive, Louisville, TN 37777 USA
TeamBest Global Companies HQ 7643 Fullerton Road, Springfield, Virginia 22153 USA
BCS tel: 865 982 0098 HQ tel: 703 451 2378 www.bestcyclotron.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA