

A Message from Krishnan Suthanthiran President & Founder, TeamBest Global Group of Companies

TeamBest Global plans to establish the technologies shown below in every empty bunker that housed a Teletherapy machine, by providing Best Total Solutions around the world!

Best Wishes,

Krishnan Suthanthiran

For more information about TeamBest Global Group of Companies, please visit www.teambest.com or contact Vidyun@teambest.com or Krish@teambest.com.

TEAMBEST GLOBAL®




www.teambest.com



Best GammaBeam & Best X-Beam External Beam Therapy Systems with IMRT/IGRT Capabilities

The Best External Beam Therapy Systems with IMRT/IGRT capabilities are an evolution of the Theratron line of treatment devices using Best GammaBeam & Best X-Beam. The advanced design provides freedom in treatment planning and can interface to all of the major record/verify systems to allow for rapid treatment parameter loading, treatment set up verification and the recording of the delivery.



NEW! Best Sub-Compact Model 200 Cyclotron with Optional Chemistry Module II (Patent Pending)

This unique low-energy, compact, fully-shielded cyclotron system is one of its kind, capable of producing F-18, C-11, N-13 and Ga-68 based radiopharmaceuticals for a variety of medical applications. The Model 200 comes standard with an F-18/FDG synthesis Best Chemistry Module I & an optional Best Chemistry Module II to produce a variety of radio tracers.

Best™ Raycell Mk1 & Mk2 Blood/Research Irradiators



Best HDR Afterloader



- Provides 20 channels for dose delivery
- Battery back-up in case of power failure and automated wire recovery
- Treatment data maintained in case of failure and treatment resumed after fix
- Quick source replacement process thereby reducing down time
- Available with Cobalt-60 or Iridium-192 sources. Also, Ytterbium-169 source available soon!

TeamBest Global Companies © 2021

* Certain products shown may not be available for sale currently.

TeamBest ASIA Private Limited email: asia@teambest.com

2/299A Plot No. 36, 5th Cross Street • V.G.P. Prabhu Nagar, Perumbakkam • Chennai – 600100 India

www.teambestasia.in www.theratronics.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA